



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J44310**

PRODUCT DESCRIPTION:	HAV P2C, Recombinant Antigen
SOURCE:	<u>E. coli</u>
PURITY:	>90% pure by SDS-PAGE (Coomassie staining)
PURIFICATION METHOD:	Inclusion bodies
BUFFER:	10 mM CBB, 0.1% SDS, 50% glycerol, pH 9.6
STORAGE TEMPERATURE:	Upon receipt store at -20°C
STABILITY:	Five years at -20°C; one month at room temperature.
SPECIFICITY:	Immunoreactive with sera of HAV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.