



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J44340**

|                      |   |
|----------------------|---|
| PRODUCT DESCRIPTION: | HAV P2C-P3A, Recombinant Antigen                          |
| SOURCE:              | <u>E. coli</u>  |
| PURITY:              | >90% pure by SDS-PAGE (Coomassie staining)                |
| PURIFICATION METHOD: | Inclusion bodies  |
| BUFFER:              | 10 mM CBB, 0.1% SDS, 50% glycerol, pH 9.6                 |
| STORAGE TEMPERATURE: | Upon receipt store at 20°C                                |
| STABILITY:           | Five years at -20°C; one month at room temperature.       |
| SPECIFICITY:         | Immunoreactive with sera of HAV-infected individuals.     |
| HAZARD/BIOHAZARD:    | There is no known biohazard associated with this product. |