



PRODUCT SPECIFICATION

CATALOG NUMBER: J44350

PRODUCT DESCRIPTION:	HAV PC3 Protein, Recombinant Antigen
SOURCE:	<u>E. coli</u>
PURITY:	>90% pure by SDS-PAGE (Coomassie staining)
PURIFICATION METHOD:	Inclusion bodies
BUFFER:	10 mM CBB, 0.1% SDS, 50% glycerol, pH 9.6
STORAGE TEMPERATURE:	Upon receipt store at 20°C
STABILITY:	Five years at -20°C; one month at room temperature.
SPECIFICITY:	Immunoreactive with sera of HAV-infected individuals.
APPLICATION:	Antigen in ELISA and Western blots, excellent antigen for detection of HAV with minimal specificity problems.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.