



PRODUCT SPECIFICATION

CATALOG NUMBER: J44370

| | |
|----------------------|---|
| PRODUCT DESCRIPTION: | HAV VP1-P2A protein, Recombinant Antigen |
| SOURCE: | <u>E. coli</u> |
| PURITY: | >90% pure by SDS-PAGE (Coomassie staining) |
| PURIFICATION METHOD: | Inclusion bodies |
| BUFFER: | 10 mM CBB, 0.1% SDS, 50% glycerol, pH 9.6 |
| STORAGE TEMPERATURE: | Upon receipt store at 20°C |
| STABILITY: | Five years at -20°C; one month at room temperature. |
| SPECIFICITY: | Immunoreactive with sera of HAV-infected individuals. |
| HAZARD/BIOHAZARD: | There is no known biohazard associated with this product. |