



PRODUCT SPECIFICATION

catalog number:	J44400
product description:	Hepatitis B Virus Core Antigen (C-Term Truncated), Recombinant
lot number:	[lot specific]
source:	Hepatitis B core protein (amino acids 1-156)
form:	Colorless to opalescent solution
purity:	100% by SDS-PAGE
buffer:	50mM Tris HCl, 500mM NaCl, 1mM EDTA, pH 7.4
preservative:	0.1% NaN ₃
molecular weight:	17 kDa by SDS-PAGE
concentration:	[lot specific] mg/mL
manufacture date:	[lot specific]
storage:	2-8°C short term (one week); <-80°C long term term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.