



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J45610**

PRODUCT DESCRIPTION:	Hsv-1 gD Antigen Recombinant aa 266-394
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	GS-4B Sepharose-Affinity purification
PURITY:	95% pure by SDS-PAGE (Coomassie staining)
BUFFER:	25 mM Tris HCl, 1mM EDTA, 50% glycerol
STORAGE TEMPERATURE:	2-8°C Short Term -20°C
STABILITY:	Fice years frozen
SPECIFICITY:	Immunoreactive with sera of Rubella HSV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.