



PRODUCT SPECIFICATION

CATALOG NUMBER: J45620

PRODUCT DESCRIPTION:	Recombinant HSV-1 Antigen protein gG
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Sepharose-derived purification
PURITY:	>95% pure by 10% SDS-PAGE (Coomassie staining)
BUFFER:	25 mM Tris HCl, 1mM EDTA, 50% glycerol
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
STABILITY:	One year frozen. One month in solution at room temperature.
SPECIFICITY:	Immunoreactive with sera of HSV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product. This product is a recombinant and is non-infectious. It is not derived from human HSV2 virus.