



PRODUCT SPECIFICATION

CATALOG NUMBER: J45630

PRODUCT DESCRIPTION:	Recombinant HSV-2 Antigen protein gD
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	GS-4B Sepharose-Affinity purification
PURITY:	95% Pure by 10% SDS-PAGE (Coomassie staining)
BUFFER:	25 mM Tris HCl, 1mM EDTA, 50% glycerol
STORAGE TEMPERATURE:	2-8°C short term. -20°C long term
SPECIFICITY:	Immunoreactive with sera of HSV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product. This product is a recombinant and is non-infectious. It is not derived from human HSV2 virus.