



PRODUCT SPECIFICATION

CATALOG NUMBER: J46000

PRODUCT DESCRIPTION:	Rubella Antigen (HPV77 strain)
SOURCE:	Tissue Culture (Vero cells)
PURIFICATION:	Purified by ultracentrifugation
PURITY:	Three major bands characteristic of Rubella antigen are detected by SDS PAGE
STORAGE TEMPERATURE:	-70°C
INFECTIVITY ASSAY:	The inactivated antigen is tested against Vero cells. The culture is followed for HA activity daily for 7 days.