



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J46020**

PRODUCT DESCRIPTION:	Rubella Antigen, 95% Pure
STRAIN:	HPV-77
PURIFICATION:	Harvested supernatant is further purified by sucrose density gradient centrifugation.
PURITY:	95%
BUFFER:	NTE buffer (Sodium chloride, TRIS, and EDTA)
PRESERVATIVE:	None
STORAGE TEMPERATURE:	-70°C long term, avoid repeated freezing and thawing
INFECTIVITY:	UV light inactivated
BIOHAZARD:	The reagent has been inactivated and tests show no infectious material. There is no known biohazard associated with this product.