



PRODUCT SPECIFICATION

CATALOG NUMBER: J46050

PRODUCT DESCRIPTION:	Rubella Virus Capsid C regions, 1-123 aa, recombinant protein
SOURCE:	<u>E. coli</u>
PURIFICATION:	Purified by proprietary chromatographic technique
PURITY:	Protein is >95% pure by 10% PAGE (coomassie staining)
BUFFER:	20 mM imidazol, 8M urea, 0.3M NaCl
STORAGE TEMPERATURE:	Shipped at ambient temperature. Upon arrival, store at -20°C.
STABILITY:	Five years frozen. One month in solution at room temperature.
SPECIFICITY:	Immunoreactive with sera of Rubella virus-infected individuals.
APPLICATION:	Antigen in ELISA and Western blots, excellent antigen for detection of Rubella Virus with minimal specificity problems.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.