



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J46060**

PRODUCT DESCRIPTION:	Rubella E1 Mosaic Antigen Recombinant gpE1 regions (157-176/374-390/213-239aa) and capsid protein
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Ion exchange chromatography
PURITY:	>95% pure by 10% PAGE (coomassie staining)
BUFFER:	20 mM imidazol, 8 M urea, 0.3 M NaCl
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
STABILITY:	Five years frozen
SPECIFICITY:	Immunoreactive with sera of Rubella virus-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.