



PRODUCT SPECIFICATION

CATALOG NUMBER: J46080

PRODUCT DESCRIPTION:	Rubella Virus E2 Antigen Recombinant
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Ion exchange chromatography
PURITY:	>95% pure by 10% PAGE (coomassie staining)
BUFFER:	20 mM imidazol, 8 M urea, 0.1 M NaCl
STORAGE TEMPERATURE:	Shipped at 2-8°C Store at -20°C
STABILITY:	Five years frozen
SPECIFICITY:	Immunoreactive with sera of Rubella virus-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.