



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J46300**

PRODUCT DESCRIPTION:	Measles Antigen
STRAIN:	Edmonston (ATCC VR-24)
FORM:	Liquid frozen at -70°C
PURIFICATION:	Antigen is prepared from the medium of infected Vero cell cultures. It is partially purified to reduce host cell components.
BUFFER:	60 mM glycine NaCl buffer, pH 9.6
PRESERVATIVE:	0.09% NaN <sub>3</sub>
STORAGE TEMPERATURE:	-70°C
APPLICATION:	Antigen performs well in solid phase immunoassays, giving high sensitivity and a low background. For in vitro use only.