



PRODUCT SPECIFICATION

CATALOG NUMBER: J46520

PRODUCT DESCRIPTION:	<i>Toxoplasma gondii</i> Antigen, crude prep, (sonicated/centrifuged)
SOURCE:	Mice (CD-1)
STARTING MATERIAL:	1x10 ⁹ tachyzoites/mL
STRAIN:	RH
BUFFER:	PBS, pH 7.4, 1 mM AEBSF, 10 KIU Aprotinin/ml
BIOBURDEN:	None
INACTIVATION:	This material is inactivated with gamma irradiated (30k rad).
INFECTIVITY CHALLENGE:	Passed
STORAGE:	≤ -70°C
HAZARD/BIOHAZARD:	This material has been inactivated . However it should be handled as potentially infectious.