



PRODUCT SPECIFICATION

catalog number:	J46550
product description:	<i>Toxoplasma gondii</i> P24 Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid
purification method:	Sepharose-affinity purification
purity:	>95% by 10% PAGE
buffer:	50mM Tris-HCl, 1.5M urea, 50% glycerol
concentration:	[lot specific] mg/mL
specificity:	This protein is immunoreactive with sera of <i>Toxoplasma gondii</i> -infected individuals.
application/notes:	This protein contains the P24 (GRA1) immunodominant region.
storage:	2-8°C (six months); -20°C (up to five years) avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.