



## PRODUCT SPECIFICATION

catalog number:	<b>J46810</b>
product description:	Zika NS1 Antigen, Recombinant
lot number:	[lot specific]
source:	Insect cells (contains His tag)
form:	Liquid
purification method:	Affinity chromatography
purity:	>85% (by SDS-PAGE)
buffer:	PBS, pH 7.4
preservative:	0.09% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL (by Bradford assay)
specificity:	This protein is reactive with IgG and IgM antibodies.
application/notes:	This protein is suitable for use in ELISA or lateral flow assays.
storage:	-20°C avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.