



PRODUCT SPECIFICATION

catalog number:	J46811
product description:	Zika NS1 Antigen, Recombinant
lot number:	[lot specific]
immunogen:	Recombinant Zika NS1 protein
source:	Insect cell expression
form:	Liquid
purification method:	Affinity purified
purity:	≥90% by SDS-PAGE
buffer:	PBS, pH 7.2
preservative:	0.02% NaN ₃
concentration:	[lot specific] mg/mL
storage:	-20°C or colder avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.