



PRODUCT SPECIFICATION

catalog number:	J46830
product description:	Zika Envelope Antigen (250-410 a.a), Recombinant, 6xHis tag at C-terminus. Strain ZikaSPH2015
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid, sterile filtered
purity:	>95% by SDS-PAGE
buffer:	PBS, 25mM arginine
concentration:	[lot specific] mg/mL
Cross-reactivity:	This protein shows minimal cross-reactivity to other Flavivirus such as Dengue and West Nile viruses
molecular weight:	~19 kDa
storage:	2-8°C short term (up to 1 week); -18°C or lower long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known biohazard associated with this product.