



PRODUCT SPECIFICATION

catalog number:	J46831
product description:	Zika Envelope Full Length Antigen (40-400 a.a), Recombinant, 6xHis tag at C-terminus
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid, sterile filtered
purity:	>95% by SDS-PAGE
buffer:	PBS, 25mM arginine
concentration:	[lot specific] mg/mL
molecular weight:	~38 kDa
storage:	2-8°C short term (up to 1 week); -18°C or lower long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known biohazard associated with this product.