



## PRODUCT SPECIFICATION

catalog number:	<b>J46835</b>
product description:	Zika Envelope Sf9 Antigen, Recombinant, 6xHis tag at C-terminus
lot number:	[lot specific]
source:	<i>Sf9, Baculovirus cells</i>
form:	Liquid, sterile filtered
purity:	>85% by SDS-PAGE
buffer:	PBS, pH 7.4 and 0.1% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
molecular weight:	~50 kDa
storage:	2-8°C short term (up to 1 week); -18°C or lower long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.