



PRODUCT SPECIFICATION

catalog number:	J46836
product description:	Zika Ectodomain Antigen, Recombinant,
lot number:	[lot specific]
source:	<i>Baculovirus Insect Cells</i>
form:	Liquid, sterile filtered
purity:	>90% by SDS-PAGE
buffer:	10mM Sodium Phosphate, pH7.2, 150mM NaCl
concentration:	[lot specific] mg/mL
molecular weight:	~55 kDa
storage:	2-8°C. Do not freeze
hazard/biohazard:	There is no known biohazard associated with this product.