



PRODUCT SPECIFICATION

catalog number:	J46840
product description:	Zika Envelope-NS1 Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i> -derived (contains C-term 6-His tag)
form:	Liquid
purity:	>90% by SDS-PAGE
buffer:	20mM phosphate, 10mM NaCl, pH 8.0
concentration:	1.0 mg/mL
specificity:	This protein is immunoreactive with sera of Zika-infected individuals.
storage:	2-8°C (shipping); -80°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.