



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>J47510</b>   |
| product description: | Streptolysin O (SLO), Recombinant Antigen                           |
| lot number:          | [lot specific]  |
| source:              | <i>E. coli</i>  |
| purity:              | ~90% by SDS Page (Commassie staining)                               |
| buffer:              | 20mM phosphate buffer, 150mM NaCl, 50% glycerol                     |
| SLO content:         | [lot specific] HU/mg (by RBC Hemolysis)                             |
| bioburden:           | < 10 CFU/mL   |
| total protein:       | [lot specific] mg/mL ( $A_{280nm}$ )                                |
| expiration date:     | [lot specific]  |
| storage:             | -20°C   |
| hazard/biohazard:    | There is no known hazard or biohazard associated with this product. |