



## PRODUCT SPECIFICATION

catalog number:	<b>J47511</b>
product description:	Streptolysin O (SLO), Recombinant Antigen (Glutaraldehyde stabilized, contains free amino group)
lot number:	[lot specific]
source:	<i>E. coli</i>
purification:	Metal affinity chromatography
buffer:	20mM Tris-HCl, 150mM NaCl, pH 8.2±0.2
preservative:	0.09% NaN <sub>3</sub>
concentration:	[lot specific] mg/ml (by OD280)
bioburden:	≤10 CFU/mL
total activity:	[lot specific] IU/mL (OD280)
specific activity:	[lot specific] IU/mL
expiration date:	[lot specific]
storage:	-70°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.