



PRODUCT SPECIFICATION

CATALOG NUMBER: J47512

PRODUCT DESCRIPTION:	Stabilized SLO, Recombinant Antigen
SOURCE:	<i>E. coli</i>
PURIFICATION:	Metal affinity chromatography
FORM:	Liquid (aseptically filtered)
BUFFER:	20 mM TRIS-HCl, 150 mM NaCl, pH 8.2 ± 0.2
PRESERVATIVE:	0.09% sodium azide
BIOBURDEN:	< 10 CFU/ml
STORAGE TEMPERATURE:	-70°C. Avoid repeated freeze/thaw cycles.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known hazard associated with this product.