



PRODUCT SPECIFICATION

catalog number:	J47512
product description:	Streptolysin O, Recombinant Antigen (Glutaraldehyde stabilized, contains free aldehyde group)
lot number:	[lot specific]
source:	<i>E. coli</i>
purification:	Metal affinity chromatography
form:	Liquid, 0.2 μ aseptically filtered
buffer:	20mM Tris, 0.15M NaCl, pH 8.2 \pm 0.2
preservative:	0.09% NaN ₃
specific activity	[lot specific] IU/mg
total activity:	[lot specific] IU/mL
bioburden:	< 10 CFU/mL
total protein:	[lot specific] mg/mL (A _{280nm})
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	-70°C. Avoid repeated freeze/thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.