



PRODUCT SPECIFICATION

| | |
|----------------------|--|
| catalog number: | J47620469 |
| product description: | HEV ORF2/ORF3 Antigen, Recombinant |
| lot number: | [lot specific] |
| source: | <i>E. coli</i> derived, protein contains immunodominant regions from ORF2 and ORF3 |
| form: | Liquid |
| purification method: | GS-4B sepharose-affinity purification |
| purity: | >95% pure by SDS-PAGE |
| buffer: | 25mM Tris-HCl, 1mM EDTA, 0.5 mM urea, 50% glycerol |
| concentration: | [lot specific] mg/mL |
| specificity: | Immunoreactive with sera of HEV-infected individuals. |
| application/notes: | This product is stable for one month at room temperature. |
| storage: | 2-8°C short term; -20°C (up to five years) |
| hazard/biohazard: | There is no known hazard or biohazard associated with this product. |