



## PRODUCT SPECIFICATION

catalog number:	<b>J47630469</b>
product description:	HEV ORF2 Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i> derived, protein contains immunodominant regions from ORF2
form:	Liquid
purification method:	GS-4B sepharose-affinity purification
purity:	>95% pure by SDS-PAGE
buffer:	25mM Tris-HCl, 1mM EDTA, 0.5 mM urea, 50% glycerol
concentration:	[lot specific] mg/mL
specificity:	Immunoreactive with sera of HEV-infected individuals.
application/notes:	This product is stable for one month at room temperature.
storage:	2-8°C short term; -20°C (up to five years)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.