



## PRODUCT SPECIFICATION

catalog number:	<b>J48040</b>
product description:	HIV-1 gag p24 Antigen, Recombinant, aa 77 – aa 436
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid
purity:	>95% by SDS-PAGE
buffer:	10 Mm PBS, Ph 7.4
molecular weight:	39 kDa
concentration:	[lot specific] mg/mL
fusions:	B-galactosidase (114 kDa) fused at the N-terminus
specificity:	Reacts strongly with human HIV positive serum
storage:	2-8°C short term; -80°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.