



PRODUCT SPECIFICATION

catalog number:	J48150
product description:	HIV-1 p24 Core Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid
purification method:	S-Sepharose
purity:	>95% by HPLC-C4 and 10% PAGE
buffer:	50 mM Tris, 8 M Urea, pH 8.5
concentration:	1 mg/mL
specificity:	The protein contains the full-length sequence of HIV-1 p24
notes/applications:	ELISA and Western blots
storage:	2-8°C short term (6 months); -20°C long term (5 years). Avoid multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.