



PRODUCT SPECIFICATION

catalog number:	J48152
product description:	HIV-2 envelope Antigen, Recombinant
product information:	The protein contains HIV-2 immunodominant regions of env gp36.
lot number:	[lot specific]
source:	<i>E. coli</i> derived 230 Amino Acids long non-fusion recombinant protein
form:	Liquid
purification method:	S-Sepharose, affinity purification
purity:	>95% by HPLC-C4 and 10% PAGE
buffer:	50 mM Tris HCL, 8 M Urea, pH 8.4,
concentration:	1.0 mg/mL
specificity:	Immunoreactive with all sera of HIV-2 infected individuals.
notes/applications:	ELISA and Western blots
storage:	2-8°C (short term one month); -20°C (up to 5 years)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.