



PRODUCT SPECIFICATION

catalog number:	J48153
product description:	HIV-2 gp36 derived Polypeptide
product information:	Chemically synthesized polypeptide containing sequences from the HIV-2 gp36 specific immune dominant domain
lot number:	[lot specific]
form:	Liquid
purification method:	S-Sepharose, affinity purification
purity:	>95% by HPLC-C18 and MASS Analysis
buffer:	H ₂ O
concentration:	1.0 mg/mL
specificity:	Immunoreactive with all sera of HIV-2 infected individuals.
notes/applications:	Rapid tests, ELISA and Western blots
storage:	2-8°C (short term one month); -20°C (up to 5 years)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.