



PRODUCT SPECIFICATION

catalog number:	J48154
product description:	HIV type O gp 41 derived peptide
product information:	Chemically synthesized polypeptide containing sequences from the HIV type O (MVP5180 and consensus sequences) gp 41.
lot number:	[lot specific]
form:	Liquid
purity:	>95% by HPLC-C18 and MASS Analysis
buffer:	H ₂ O
concentration:	1.0 mg/mL
specificity:	Immunoreactive with all sera of HIV type O infected individuals.
notes/applications:	Rapid tests, ELISA and Western blots
storage:	2-8°C (up to one year)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.