



PRODUCT SPECIFICATION

catalog number:	J48157501
product description:	HIV gp-24 Antigen, Purified
lot number:	[lot specific]
source:	<i>In vitro</i> purified proteins from <i>E. Coli</i>
form:	liquid
buffer:	50mM Borax Buffer, pH 9.0
preservative:	0.02% NaN ₃
purity:	>90% SDS-PAGE
concentration:	[lot specific] mg/mL, absorbance at 280 and 260 nm
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.