



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J48270**

PRODUCT DESCRIPTION:	PR3 Antigen, High Pure
SOURCE:	Human Neutrophils
PURITY:	>98% by SDS-PAGE
FORM:	Liquid frozen
BUFFER:	50 mM sodium acetate, 0.1 M NaCl, 0.02% lubrol, pH 4.5 ± 0.2
TOTAL PROTEIN:	0.211 mg/ml (Bradford method)
BIOBURDEN:	< 10 CFU/ml
STORAGE TEMPERATURE:	-20°C. Avoid repeated freeze/thaw cycles.
HAZARD/BIOHAZARD:	All donor units have been tested and found negative for HBsAg, anti-HCV and HIV I & II antibodies by FDA approved methods; found negative for viral DNA/RNA for HBV, HIV <sub>1</sub> and HCV by PCR. However, this material should still be handled as potentially hazardous.