



## PRODUCT SPECIFICATION

catalog number:	<b>J49080</b>
product description:	HIV-II gp160 Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i> derived, N-terminal $\beta$ -gal fusion
form:	Liquid
purity:	>95% by SDS-PAGE
buffer:	20mM Tris-HCl, 8M Urea, 10mM $\beta$ ME, pH 8.0
concentration:	[lot specific] mg/mL
notes/applications:	This antigen reacts strongly with human HIV positive sera.
storage:	2-8°C (< 3 months); -80°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.