



PRODUCT SPECIFICATION

CATALOG NUMBER: J49400

PRODUCT DESCRIPTION:	Cystatin C Antigen, Part Pure
SOURCE:	Human urine
PURITY:	10-50%
FORM:	Homogeneous liquid (aseptically filtered)
BUFFER:	10 mM Tris, 10 mM Benzamidine, 20% Glycerol, pH 8.5 ± 0.2
PRESERVATIVE:	0.09% sodium azide
CYSTATIN C CONTENT:	1-2 mg/ml (by Roche Tina quant)
BIOBURDEN:	< 10 CFU/ml
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	All donor units have been tested and found negative for HbsAg, anti-HCV and HIV I & II antibodies by FDA approved methods; found negative for viral DNA/RNA for HBV, HIV ₁ and HCV by PCR. This material contains sodium azide as a preservative, take appropriate precaution.