



## PRODUCT SPECIFICATION

catalog number:	<b>J50001</b>
product description:	Human SAA (Ro-52) Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid
purity:	>90% (by SDS-PAGE)
buffer:	20mM phosphate buffer, 500mM NaCl, 3M Urea, pH 6
concentration:	1.0 mg/mL
storage:	-70°C or below avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.