



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J50010**

PRODUCT DESCRIPTION:	Recombinant human Ro-60 (SSA)
SOURCE:	Expressed in <u>E. coli</u>
BUFFER:	6 M urea, 500 mM NaCl, 10 mM Tris, 500 mM Imidazole (pH 6.0)
SDS PAGE:	Major 60 kDa band
WESTERN BLOT:	60 kD band strongly reactive with human anti-Ro-SSA antisera
STORAGE CONDITIONS:	-20 to -70°C (avoid repeated freeze thaw cycles)
COMMENTS:	To coat ELISA plates we suggest at least a 100-fold dilution in 50 mM carbonate buffer (pH 9.6). Ensure that the sample is mixed and used immediately.