



PRODUCT SPECIFICATION

CATALOG NUMBER: J50020

PRODUCT DESCRIPTION:	C-terminal truncation of the human La/SSB autoantigen with an N-terminal TnC tag
SOURCE:	Expressed in <u>E. coli</u>
BUFFER:	10 mM Sodium phosphate, 150 mM NaCl, 0.05% NaN ₃ (pH 7.5)
SDS PAGE:	Major 47 kDa band
STORAGE:	-70°C (avoid repeated freeze thaw cycles)
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.