



PRODUCT SPECIFICATION

CATALOG NUMBER: J50030

PRODUCT DESCRIPTION:	Recombinant Scl-70 Antigen
SOURCE:	Expressed in Baculovirus-infected Sf9 cells
PURIFICATION:	Purified by Affinity Chromatography
BUFFER:	6 M HEPES/KOH, 400 mM NaCl, pH 7.5, 20% Glycerol
SDS PAGE:	Band at 109 kD
STORAGE:	-20 to -70°C (avoid repeated freeze thaw cycles)