



PRODUCT SPECIFICATION

CATALOG NUMBER: J50040

| | |
|----------------------|---|
| PRODUCT DESCRIPTION: | Human Jo-1 with an N-terminal TnC tag |
| SOURCE: | Expressed in <u>E. coli</u> |
| BUFFER: | 6 M urea, 500 mM NaCl, 10 mM Tris, pH 8.0 |
| SDS PAGE: | Major ~55 kDa band, > 95% pure |
| STORAGE: | -20°C to -70°C (avoid repeated freeze thaw cycles) |
| HAZARD/BIOHAZARD: | There is no known biohazard associated with this product. |