



PRODUCT SPECIFICATION

CATALOG NUMBER: J58150

PRODUCT DESCRIPTION:	Soluble Transferrin Receptor Antigen, Part Pure
SOURCE:	Human plasma
BUFFER:	25 mM Tris, 0.15 M NaCl, pH 8.0 ± 0.2
PRESERVATIVE:	0.09% sodium azide
BIO-BURDEN:	< 10 CFU/ml
STORAGE TEMPERATURE:	-20°C
HAZARD/BIOHAZARD:	Donor materials were tested and found negative for HBsAg, anti-HCV and anti-HIV 1 & 2 by FDA approved methods. Donor material was also tested and found negative for viral DNA/RNA for HIV 1, HCV, and HBV by PCR. However, this material should still be handled as potentially hazardous. This material contains sodium azide as a preservative, take appropriate precaution.