



## PRODUCT SPECIFICATION

catalog number:	<b>J58311</b>
product description:	Human Serum Amyloid A (SAA), Recombinant Antigen
lot number:	[lot specific]
source:	<i>E. coli</i>
form:HAZ	Lyophilized white powder, sterile filtered
purity:	≥98% (by SDS-PAGE and HPLC)
buffer:	20mM Tris-HCl, pH 9.0, 150mM NaCl
molecular weight:	~11.7 kDa
concentration;	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C (shipping); and 1 month under sterile conditions after reconstitution  -20°C to -70°C (12 months from date of receipt as supplied)  -20°C to -70°C (3 months under sterile conditions after reconstitution)
hazard/biohazard:	There is no known hazard or biohazard associated with this product.