



PRODUCT SPECIFICATION

CATALOG NUMBER: J63000

PRODUCT DESCRIPTION:	PSA Antigen
SOURCE:	Human Seminal Fluid
FORM:	Liquid in a buffer system
PURITY:	>96% by SDS-PAGE
BUFFER:	0.05M phosphate, 0.15M NaCl, pH 7.5, filter sterilized through a 0.2 uM membrane.
PRESERVATIVES:	0.1% NaN ₃
STORAGE TEMPERATURE:	2-8°C
BIOHAZARD:	Donor level tested negative for HBsAg, HIV I&2 antibody, and HCV antibody by FDA approved methods.