



PRODUCT SPECIFICATION

CATALOG NUMBER: J63020

PRODUCT DESCRIPTION:	Human Alpha-2 Macroglobulin
SOURCE:	Human serum
PURITY:	>96%
SDS-PAGE:	One band corresponding to alpha 2M molecular weight of 725kD
FORM:	Lyophilized
BUFFER:	0.02M TRIS, 0.15M NaCl and 0.1M sucrose
STORAGE CONDITION:	2-8° C
NOTE:	Starting material tested negative for HIV 1&2 antibodies, HCV antibody, and HBsAg by FDA approved methods.