



## PRODUCT SPECIFICATION

catalog number:	<b>J63020</b>
product description:	Alpha-2-Macroglobulin Antigen
lot number:	[lot specific]
source:	Human plasma
form:	Liquid
purity:	>95% by SDS-PAGE
buffer:	Frozen in 14.6ml 30 mM Na phosphate, pH 7.0 and 1mg glycine per 1mg A2M.
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	≤-20°C
hazard/biohazard:	Donor material has been tested and found non-reactive for HBsAg, anti-HCV, anti-HBc and negative for anti-HIV 1&2 by FDA approved methods. However, this material should still be handled with caution as potentially hazardous.