



PRODUCT SPECIFICATION

catalog number:	J63710
product description:	Malaria HRP-2 Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid
purity:	≥95% by SDS-PAGE
buffer:	50 mM Tris, 100 mM Sodium Chloride, 0.2M Sucrose, pH 8.0
concentration:	[lot specific] mg/mL
pairing:	This antibody forms a match pair with A63820314P, A63830314P and A63840314P.
storage:	-20°C avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.