



PRODUCT SPECIFICATION

CATALOG NUMBER: J64020

PRODUCT DESCRIPTION:	Beta-2-Microglobulin Antigen, High Pure
SOURCE:	Human Urine
FORM:	Liquid
BUFFER:	40 mM PBS, 1 mM EDTA, 0.15 M NaCl, pH 7.2
PRESERVATIVE:	0.05% sodium azide
PURITY:	>98% by SDS-PAGE, single band ~12kD
STORAGE:	2-8°C short term; -20°C long term.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. All donor units were tested and found non-reactive for HBsAg, anti-HCV and anti-HIV 1&2 by BioRad and OCD. Donor units tested and found negative for viral HBV, HCV and HIV 1 by PCR.